

Spectrum IQ Infusion System

Integrating the Healthcare Enterprise (IHE)

Integration Statement

The IHE Integration Statement below describes the intended conformance of the **Spectrum IQ** Infusion System with the IHE Technical Framework.

| Vendor | Product Name | Date |
|---|--|--|
| Baxter | Spectrum IQ Gateway | January 2018 |
| <i>This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:</i> | | |
| Profiles Implemented | Actor Implemented | Options Implemented |
| Device Enterprise Communications (DEC) | Device Observation Reporter (DOR) | PCD-01 (Communicated PCD Data) |
| Infusion Pump Event Communications (IPEC) | IPEC - DOR (Device Observation Reporter) | PCD-10 (Communicate Infusion Pump Event) |
| Point-of-Care Infusion Verification (PIV) | Infusion Order Consumer | PCD-03 (Communicate Infusion Order) |
| Alert Communication Management (ACM) | Alert Reporter (AR) | PCD-04 (Report Alert) |
| Consistent Time (C1) | Time Client | ITI TF-1 |
| Medical Equipment Management Location Services (MEMLS) | Location Observation Consumer (LOC) | PCD-16 (Report Location Observation) |
| Links to general information on IHE | | |
| In North America: www.ihe.net | | Profiles: www.ihe.net/Patient_Care_Devices/ |

Rx Only. For safe and proper use of this device, refer to the appropriate operator's manuals.

www.baxter.com
www.SpectrumIQ.com

Baxter Healthcare Corporation
 One Baxter Parkway / Deerfield IL 60015
 Baxter and Spectrum IQ are trademarks of Baxter International Inc.